



Copal Cavity Varnish with Fluoride

Protective varnish based on Copal gum for insulating exposed dentine and providing a thermal barrier.

Indications

- Ideal cavity liner and dentinal tubule sealer
- Assures perfect protection to exposed dentine
- Intermediary varnish prior to insertion of amalgams, cements or composites

Product Benefits

- Insulates exposed dentine
- Reduces penetration of fluids into the tooth
- Fluoride provides a thermal barrier for deep restoration
- Fluoride alleviates thermal sensitivity
- Needs no thinners
- Compatible with most amalgams, cements and composites
- High fluidity for easy application