



Non-Eugenol Temporary Luting Cement

OraTemp NE is a non-eugenol temporary cement that permits placement and retention under the function of permanent crowns, bridges, splints and their temporary counterparts. OraTemp NE is insoluble in oral fluids and retains a tight marginal seal to prevent leakage of these fluids under the restorations.

Indications

- Cementation of crowns, bridges and splints
- Retentions of temporary crowns and bridges
- Trial cementation of permanent restorations prior to final insertion

Product Benefits

- Excellent physical properties
- Easy to mix – easy to clean
- Low solubility in oral fluids
- Increased patient comfort

Product Packaaina